

dNerva™

Targeted Lung Denervation

A one-time bronchoscopic procedure designed to reduce the frequency and severity of COPD exacerbations



Nuvaira® Lung Denervation System is an investigational device in the United States and has CE mark regulatory approval.

NUVAIRA®

Understanding the “Lung Attack”

Chronic obstructive pulmonary disease (COPD) is progressive and leaves patients experiencing symptoms of breathlessness, coughing, sputum production, chest tightness and wheezing driven by the natural release of acetylcholine.

The brain signals the lung, via the vagus nerve, to release acetylcholine in response to toxins, such as tobacco smoke or environment pollution. This generates a repeating cycle of signaling along the nerve pathways between the lung and the brain.

In some COPD patients, this constant cycle of nerve signaling and chronic airway inflammation, creates a state of heightened neural sensitivity called airway hyperresponsiveness. This puts the pulmonary system into a state of permanent high alert which results in an over-exaggerated reaction to normal environmental stimuli. These patients may then suffer from frequent exacerbations, or lung attacks, which medications cannot adequately prevent.³

Treatment goals focus on minimizing the impact of lung attacks¹⁹, which can cause:

-  Accelerated decline in lung function^{12,13,14}
-  Increased risk of acute cardiovascular events¹⁶
-  Reduced activities of daily living (ADL) and quality of life¹⁵
-  Increased risk of future hospital admissions⁶
-  Increased mortality risk^{17,18}

Lung attacks not only reduce patient quality of life; one severe exacerbation can leave a patient at 14 times greater risk of death within one year.^{17,18}

COPD is the biggest driver of health care costs related to hospitalization and responsible for \$72 billion in annual direct health care expense.^{6,7,8,9}

Current Therapies Leave Significant Unmet Needs

Prescribed medications are an integral component of effective disease management, but patients adhering to optimal medical therapy continue to exacerbate in large-scale clinical trials targeting reduction in COPD exacerbation as the primary goal.^{10,11}

Disrupting Pulmonary Nerve Input

The goal of **Targeted Lung Denervation (TLD)** is to disrupt parasympathetic pulmonary nerve input to the lung to reduce the clinical consequences of neural hyperactivity^{1,2}. Performed during a single bronchoscopic procedure, instead of via inhaled medications, this denervation of the lung permanently disrupts the nerve signaling to and from the lung, which allows airway smooth muscle to relax and mucous secretion to decrease. This moves the pulmonary system to a state of “low alert” to improve symptom tolerance, decrease symptom burden and decrease exacerbations.

This one-time surgical procedure involves passing a specialized catheter through a flexible bronchoscope to complete a full circumferential ablation in the

mainstem bronchus of each lung. The bronchoscope is inserted into the mainstem bronchi, and a radio-opaque esophageal balloon is temporarily inserted to facilitate visualization and avoidance of treatment near the esophagus. Performed under general anesthesia, the procedure typically takes 60–90 minutes.

The ablation permanently disrupts pulmonary nerve input to the lung, reducing acetylcholine release. With less acetylcholine, the airway smooth muscle relaxes, and mucous secretion decreases, moving the pulmonary system to a state of “low alert” to improve symptom tolerance, decrease symptom burden and decrease exacerbations.

The Procedure

Step 1: Position

The dNerva catheter threads through the working channel of the bronchoscope and is positioned in the distal aspect of the mainstem bronchi.

Step 2: Inflate

The conduit and balloon are inflated by circulating coolant through the catheter which brings the electrode in contact with the airway wall.

Step 3: Confirm


Visual confirmation of good contact of the electrode with the airway wall and adequate distance from the esophagus (confirmed under fluoroscopy).

Step 4: Active


Radio-frequency energy activation. After completion of the activation, the physician will rotate the catheter 90 degrees and then repeat. Each mainstream bronchus is treated four times to complete circumferential ablation and thus deliver a full treatment.

Post-Procedure

The patient can expect to go home the same day as the procedure. The patient will be provided with a card explaining their participation in a clinical trial and are advised to present that should they become ill for any reason and seek medical care.



Nuvaira Clinical Development Program




1. Valipour, Chron Obstruct Pulm Dis 2018. 2. Valipour, Respiration 2019. 3. Pison, Respir Res 2021. 4. Slebos, Shah, AJRCCM 2019. 5. Valipour, Int J COPD 2020.

AIRFLOW 2 Trial Results to Date

AIRFLOW 2 was a multicenter, 1:1 randomized, sham-controlled, double-blinded, prospective study. Unblinded after 12.5 months, sham patients were invited to cross-over to TLD treatment in a separate protocol. Intended to evaluate the safety and impact of TLD in moderate-to-severe symptomatic patients on optimal medical therapy, it was offered at 16 study sites in 5 countries.

Primary Endpoint




TLD Impact on Two Year Outcomes²⁰



Observed Clinical Stability after Three Years^{21,22}



AIRFLOW 2 Results Provided Solid Basis for AIRFLOW 3



The Nuvaira Lung Denervation System is the first pulmonary product to successfully apply to France's National Authority for Health (HAS) Forfait Innovation program. The program grants temporary funding to evaluate CE marked devices that are deemed to be truly innovative in a rigorous clinical trial to confirm a significant clinical benefit, filling an unmet medical need and leading to a simplified market access pathway.

The Nuvaira Lung Denervation System has been designated as a Breakthrough Device by the US Food and Drug Administration (FDA). This program creates an expedited pathway for prioritized FDA review of devices that have potential to provide more effective treatment or diagnosis of life-threatening or irreversibly debilitating diseases or conditions. In March 2021, the FDA's review of safety data on the first 50 patients enrolled in Nuvaira's AIRFLOW 3 pivotal trial resulted in full IDE approval for completion of the 400-patient study, and in April, the Centers for Medicare & Medicaid Services (CMS) approved the AIRFLOW 3 trial for the purposes of Medicare coverage (42 CFR 405 Subpart B).

AIRFLOW 3 Is Now Enrolling

The AIRFLOW 3 trial will be comprised of up to 400 patients across sites in the United States and Europe with the primary endpoint being a comparison of the probability of subjects having a moderate or severe exacerbation between the treatment arm and the sham-control group at one year post procedure.

INCLUSION CRITERIA

40+ years old

BMI 18–35

25–80% FEV₁ predicted and FEV₁ / FVC < 70% post-bronchodilator

2+ moderate (OCS and/or antibiotics treatment) or 1+ severe (ER or hospitalization) COPD exacerbations in the past year

CAT ≥ 10

Non-smoking / Non-vaping for 2+ months

Candidate for bronchoscopy and general anesthesia

Current flu vaccine (or willingness to be vaccinated when available)

Resting SpO₂ ≥ 89% on room air

EXCLUSION CRITERIA

Ongoing GI-related symptoms (severe uncontrolled GERD or severe dysphagia, gastroparesis, or recent major GI surgery)

Diagnosis of active non-COPD lung disease

Clinically relevant bronchiectasis

Current diagnosis of active asthma

MI < 6 months, EKG evidence of life-threatening arrhythmias or acute ischemia

Pre-existing documented LVEF < 40%, stage C or D ACC/AHA or Class III or IV (NYHA) CHF

Diagnosed pulmonary hypertension with clinical evidence of CV function impairment (peripheral edema, sustained mPAP > 25 mmHg at rest (right hearth cath), or estimated RVSP > 50 mmHg echo)

Radiation or chemo treatment within 1 year of consent

> 10 mg daily prednisone or equivalent

Pulmonary Rehab: Not required for this trial. If patient is currently in rehab, wait 3 months after completion before consenting.

Endobronchial Valves: Patients may be eligible if EBV implant > 1 year. Or EBV explant > 3 months and implant > 1 year.

Opioids: Low dose use allowed for valid medical condition.

Trial Locations Include

TREATMENT LOCATIONS IN THE UNITED STATES

FirstHealth of the Carolinas
Pinehurst, NC

University of Michigan Health
Ann Arbor, MI

Washington University Medical Center
St. Louis, MO

Medical University of South Carolina
Charleston, SC

Beth Israel Lahey Health
Burlington, MA

Suburban Lung Associates
Elk Grove Village, IL

The Ohio State University Wexner Medical Center
Columbus, OH

Temple Health Lung Center
Philadelphia, PA

Ascension St. Vincent's
Jacksonville, FL

University of Louisville Health
Louisville, KY

Houston Methodist Hospital
Houston, TX

University of Pittsburgh Medical Center
Pittsburgh, PA

Spectrum Health
Grand Rapids, MI

Duke University Medical Center
Durham, NC

Harbor-UCLA Medical Center
Torrance, CA

UC Davis Health
Sacramento, CA

University of Alabama Lung Health Center
Birmingham, Alabama

St. David's Medical Center
Georgetown, TX

The University of Chicago Medicine
Chicago, IL

Honor Health
Phoenix, AZ



TREATMENT LOCATIONS IN EUROPE

Klinik Floridsdorf
Vienna, Austria

Thoraxklinik Universitätsklinikum Heidelberg
Heidelberg, Germany

CHU de Grenoble
Grenoble, France

CHU de Reims
Reims, France

Nouvel Hôpital Civil
Strasbourg, France

Claude-Bernard - Hôpital Bichat
Paris, France

Hôpital Larrey
Toulouse, France

Hôpital Arnaud de Villeneuve
Montpellier, France

Hôpital Pasteur
Nice, France

Hôpital de la Cavale Blanche
Brest, France

Royal Brompton Hospital
London, UK

University Medical Center Groningen
Groningen, Netherlands

Amsterdam UMC
Amsterdam, Netherlands

References

1. Hummel et al 2019 Physiologic and histopathologic effects of targeted lung denervation in an animal model *J Appl Physiol* 126:67-76.
2. Valipour et al 2020 Two-year outcomes for the double-blind, randomized, sham-controlled study of targeted lung denervation in patients with moderate to severe COPD: AIRFLOW 2 *Int J Chron Obstruct Pulmon* 15:2807-2816.
3. Pavord et al 2017 Mepolizumab for eosinophilic chronic obstructive pulmonary disease *N Engl J Med* 377(17): 1613-1629.
4. Mayse et al 2020 Targeted lung denervation in sheep: durability of denervation and long-term histologic effects on bronchial wall and peribronchial structures *Respir Res* 21:117.
5. Zhang et al 2018 A systematic review of how patients value COPD outcomes *Eur Respir J* 19;52(1).
6. Ihenacho et al 2020 Economic burden of chronic obstructive pulmonary disease (COPD): A systematic literature review *Int J Chron Obstruct Pulmon Dis* 26;15: 439-460.
7. Dalal et al 2011 Costs of COPD exacerbations in the emergency department and inpatient setting *Respir Med* 105(3):454-60.
8. Dalal et al 2015 Impact of COPD exacerbation frequency on costs for a managed care population *J Manag Care Spec Pharm* 21(7): 575-83.
9. Patel et al 2018 COPD affects worker productivity and health care costs *Int J COPD* 13:2301-2311.
10. Lipson et al 2018 Once-daily single-inhaler versus dual therapy in patients with COPD *N Engl J Med* 378(18): 1671-1680.
11. Rabe et al 2020 Triple inhaled therapy and two glucocorticoid doses in moderate-to-very severe COPD *N Engl J Med* 383(1):35-48.
12. Donaldson et al 2002 Relationship between exacerbation frequency and lung function decline in chronic obstructive pulmonary disease *Thorax* 57(10):847-52.
13. Halpin et al 2017 Effect of a single exacerbation on decline in lung function in COPD *Respir Med* 128: 85-91.
14. Dransfield et al 2017 Acute exacerbations and lung function loss in smokers with and without Chronic Obstructive Pulmonary Disease *Am J Respir Crit Care Med* 195(3): 324-33.
15. Seemungal et al 1998 Effect of exacerbation on quality of life with chronic obstructive pulmonary disease *Am J Respir Crit Care Med* 157:1418-22.
16. Westerik et al 2017 Associations between chronic comorbidity and exacerbation risk in primary care patients with COPD *Respir Res* 18(1):31.
17. Soler-Cataluna et al 2005 Several acute exacerbations and mortality in patients with chronic obstructive pulmonary disease *Thorax* 60(11):925-31.
18. Rothnie et al 2018 Natural history of chronic obstructive pulmonary disease exacerbations in a general practice-based population with chronic obstructive pulmonary disease *Am J Respir Crit Care Med* 198(4):464-471.
19. Global Strategy for the Diagnosis, Management and Prevention of Chronic Obstructive Pulmonary Disease 2021 Report.
20. Hurst et al 2020 Understanding the impact of chronic obstructive pulmonary disease exacerbations on patient health and quality of life *Eur J Intern Med* 73:1-6.
21. Pison et al 2021 Safety of denervation following targeted lung denervation therapy for COPD AIRFLOW 1 3-year outcomes *Respir Res* 22(1): 62.
22. Hartman JE, Conway F, Degano B, Augustijn SWS, Herth FJF, Mayr A, Chacaroun S, Tonkin J, Valipour A, Slebos DJ. Rate of lung function decline slows in the 3 years after targeted lung denervation in COPD. *Respir Med* 2021; 188: 106604.

Common potential risks associated with the TLD procedure include, but are not limited to, worsening of COPD symptoms (shortness of breath, increased cough, COPD exacerbations), coughing up blood tinged mucous, difficulty swallowing, chest pain, upset stomach or feeling of fullness, trouble processing food, and fever.

Nuvaira® Lung Denervation System is an investigational device in the United States and has CE mark regulatory approval.



Nuvaira, Inc.
6500 Wedgwood Road North, Suite 100
Maple Grove, MN 55311
Tel: 763.450.2800
www.nuvaira.com